Instructions for Use

Instructions for Use - GB
Gebrauchsanweisung - DE
Instructions d'utilisation - FR
Bruksanvisning - SE
Istruzioni per l'Uso - IT
Instrucções para a utilização - PT
Brugsanvisning - DK
Käyttöohje - FI
Manual de instrucciones - ES
Οδηγίες χρήσης - GR
Návod k použití - CZ
Instrukcja obsługi - PL
Gebruiksaanwijzing - NL
Bruksanvisning - NO
Kasutusjuhend - EE
Használati útmutató - HU
Novadila za uporabo - SL

Pediatric Audiometer PA5

***Valid from serial 01
80691304 – 10/2010***

Instructions for Use - GB

Pediatric Audiometer PA5

Intended Use

The PA5 handheld paediatric screening audiometer is designed to be a device for screening for hearing loss primarily amongst children. Output and specificity of this type of device are based on the test characteristics defined by the user, and may vary depending on environmental and operating conditions. The screening for hearing loss using this kind of audiometer depends on the interaction with the patient. However, for children not responding well, various test possibilities allow the tester of having at least some evaluative result. Thus, a "normal hearing" result should not allow for ignoring other contra indications in this case. A full audiologic evaluation should be administered if concerns about hearing sensitivity persist.

Precautions

| | |
|--|--|
|  WARNING | WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury. |
|  CAUTION | CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury. |
|  NOTICE | NOTICE is used to address practices not related to personal injury. |

NOTICE


1. Be sure to use only stimulation intensities, which will be acceptable for the patient.
2. The transducers (headphones, bone conductor, etc.) supplied with the instrument are calibrated to this instrument - exchange of transducers require a recalibration.
3. It is recommended that parts which are in direct contact with the patient (e.g. earphone cushions) are subjected to standard disinfecting procedure between patients. This includes physically cleaning and use of a recognised disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.
4. Please note that the CE Marking is only legal if this Instruction is translated into the national language of the user no later than at the delivery to him, if the national legislation demands a text in the national language according to MDD article 4.4.
5. Although the instrument fulfils the relevant EMC requirements precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.
6. Disposal of batteries must be made according to national regulations.

If this apparatus is connected to one or more other devices with medical CE marking, to make up a system or pack, the CE marking is only valid also for the combination if the supplier has issued a declaration stating that the requirements in the Medical Device Directive article 12 are fulfilled for the combination

Operation Instruction

The instructions included in this manual describe the general functions of the instrument. A more detailed discussion is provided in the Operation Manual in English.

Description of Control Panel

Indication of stimulus mode:

Indication LEDs informing the user of the present stimulus mode: Tone, NB or WN.

Light:

Light button to control the three red LEDs, which are arranged in a triangle above the speaker in order to condition the orientation reflex.

Stimulus:

Stimulus button to present the selected stimulus Tone, NB or WN.

Select:

Select button to select between the three different stimuli Tone, NB or WN.


Intensity dB HL:

With the intensity control switch intensities between 20 and 80 dB HL can be selected in steps of 10 dB when the distance between the ear and the loudspeaker of the PA5 is 50 cm or the PA5 is switched off by leaving the intensity control switch in the "Off" position. When the PA5 is not activated for two minutes it will switch off automatically.

Frequency kHz:

With the frequency control switch it is possible to select between the following frequencies: 0.5, 1, 2, 3, and 4 kHz.

Description of Stimulus Panel

Loudspeaker:

The loudspeaker is to be found underneath the black grid. When used on a patient the grid should be positioned in a distance of 50 cm (20 inches) from the ear in order to obtain the intensities printed on the Control Panel.

LEDs:

Three LEDs arranged in a triangle for conditioning of the orientation reflex.

Headphone Connector:


Connector for single Headphone TDH39S (optional). When the headphone is connected PA5 will automatically go to Pure Tone stimulation and correct calibration for Pure Tone Audiometry will be applied.

Battery Description

Replacing Batteries:

To replace old batteries unscrew the small black lid in the narrow end of PA5 and the batteries can be taken out.

Replace with three new AA batteries. When inserting the new batteries make sure that they are inserted correctly according to the small drawing in the bottom part of the control panel.

PA5 contains 3 batteries, size LR6, AA or Mignon.

Approximate Battery Lifetime:

The battery lifetime using Alkaline battery type:

| | |
|--|-----------|
| With the instrument switched off: | 12 months |
| With 80 dB tone switched on: | 10 hours |
| With 80 dB tone and light switched on: | 4 hours |

Battery Level indication:

When the batteries need to be replaced the LED indication for the present used stimulus will gradually reduce in light intensity and finally switch off.

Note: Always remove the batteries when the instrument is left unused for a longer period.

Technical Specifications

The technical specifications provided here cover the general aspects of the instrument, while more specific details can be found in the Operation Manual in English.

Standards:

IEC 60645-1, Type 5
Sound Pressure Level of the loudspeaker: ISO 389-7
Sound Pressure Level of the headphone: ISO 389-1

Medical CE-mark:

The CE-mark indicates that Interacoustics A/S meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. Approval of the quality system is made by TÜV – identification no. 0123.

Power:

Batteries : 3 x AA, LR6 or Mignon.
Alkaline or rechargeable (NiMH or NiCa).

Frequencies:

500, 1000, 2000, 3000, 4000 Hz.

Stimuli:

Warble Tone, NB, WN.

Intensities:

Distance 50 cm : 30 – 80 dB HL in 10 dB steps, Warble Tone and WN
20 – 70 dB in NB
Distance 16 cm : increases the indicated intensity by 10 dB.

Warble Frequency:

5 Hz, \pm 5%.

Stimulation using TDH39S:

Pure Tone: 500, 1000, 2000, 3000 and 4000 Hz.
Intensities of 30 – 80 dB (independent calibration applied when TDH39 is connected).

Sound Source:

Built in loudspeaker or headphone TDH39S (independent calibration registers).

Light Stimulation:

3 LED's arranged in a triangle, flash speed 5 Hz.

Tone and Light Stimulation:

Silent touch switch.

Approximate Battery Lifetime:

The battery lifetime using Alkaline battery type:
With the instrument switched off: 12 months
With 80 dB tone switched on: 10 hours
With 80 dB tone and light switched on: 4 hours

Dimensions:

L x W x H : 25 x 7 x 5 cm / 9.8 x 2.7 x 2 inches
Weight : 0.36 kg incl. batteries

Detachable Parts

- 3 Batteries
- Carrying Bag
- Operation Manual
- CE manual

Concerning repair

1. Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:

assembly operations, extensions, readjustments, modifications or repairs are carried out by authorized persons,

a 1 year service interval is maintained

the electrical installation of the relevant room complies with the appropriate requirements, and

the equipment is used by authorized personnel in accordance with the documentation supplied by Interacoustics.

2. It is important that the customer (agent) fills out the RETURN REPORT every time a problem arises and sends it to Interacoustics, Drejervænget 8, DK-5610 Assens, Denmark. This should also be done every time an instrument is returned to Interacoustics. (This of course also applies in the unthinkable worst case of death or serious deterioration to patient or user).
3. When instrument fuses need renewal, the correct type as stated on the instrument shall be used.

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