




Optical Flow

Module tutorial



Table of contents

1.	GENERAL	3
1.1.	Description	3
1.2.	Required accessories.....	3
1.3.	Patient setup	3
2.	Settings	3
2.1.	Session settings	4
2.1.1.	Presets	4
2.1.2.	Environment	4
2.1.3.	Speed.....	6
2.1.4.	Rotation	7
2.1.5.	Random direction.....	7
3.	Module.....	8
3.1.	Session	8
3.2.	Shortcuts	9
3.3.	Results.....	9
3.3.1.	Summarized results	11
3.3.1.	Graphs and report	11



1. GENERAL

1.1. Description

Optical Flow is a software used for balance disorder rehabilitation using linear or curved visual scrolling.

The patient is placed in a tunnel that the therapist can move forward, back or even turn.

1.2. Required accessories

Headset only **OR** Headset and StaticVR force platforms or MotionVR dynamic posturography platform.

1.3. Patient setup

These are simple, limited numbers of tips drawn from the literature.

Except in special cases, prefer the standing position. The use of a foam block can direct the re-education towards more "posture/proprioception" type work, mainly with the "Advanced" environment by varying the turn tilts.




RECOMMENDATION

Stay in contact with the patient who is at risk of falling during vertical stimulations, even if they are standing on a dynamic posturography platform.

2. Settings



2.1. Session settings



2.1.1. Presets

Several different protocols are available. Each protocol corresponds to specific parameters.



RECOMMENDATION

Start with a "test session" to measure patients' tolerance to the stimulation and the proposed VR environment.

For each proposed protocol, the different options are already set. Just validate the selected protocol to launch the application.

2.1.2. Environment



Environment

Environment choice

Simple Advanced

☒ Furrows

Opacity filter 0.0

☒ Enable fixed point

Spot size 30 °

Environment choice:

- **Simple:** radial optical flow, linear, without turn possibilities. 3D elements (asteroids) can be activated/deactivated using the "Q" key (See Shortcuts).
- **Advanced:** used to orientate the tunnel with bends to strengthen balance, proprioception, and work on the patient's posture.

Furrows:

Enabling this setting is possible for the "**advanced**" tunnel.

Used to have a less smooth texture of the scenery.

Opacity filter:

Allows the stimulation to appear progressively for the most sensitive patients.

Varies the transparency of the environment.

Value: 0,0 to 1,0.

Enable fixed point and point size:

Allows a fixed point to appear in the patient's field of vision.

The point's size can be adjusted using the slider.

Value: from 5 to 90°.



2.1.3. Speed

A screenshot of a software interface titled 'Speed'. It contains three sliders and a checkbox. The first slider is labeled 'Speed' and has a value of '15 m/s'. The second slider is labeled 'Forward flow duration' and has a value of '10 s'. The third slider is labeled 'Backward flow duration' and has a value of '8 s'. A checkbox labeled 'Automatic inversion' is checked.

Speed

Speed 15 m/s

☒ Automatic inversion

Forward flow duration 10 s

Backward flow duration 8 s

Speed:

Scenery backward or forward scrolling speed.

Value: -50 to 50 m/s.

Automatic inversion:

Used to program changes in the scrolling direction by indicating a **forward** and **backward flow** duration.

Forward flow duration:

Value: 10 to 120 s.

Backward flow duration:

Value: 5 to 30 s.

**ATTENTION**

The inversion causes significant anteroposterior postural reactions, and a high risk of falling.



2.1.4. Rotation

Rotation

Rotation

0 °/s

☒ Automatic inversion

Rotation duration

90 s

Rotation:

Backward or forward scenery rotation. By convention, anticlockwise rotations are negative and clockwise rotations are positive.

Automatic inversion allows the direction of rotation to be programmed by indicating a duration.

Value: -30 to 30°/s.


Automatic inversion:

The automatic reversing function allows you to program changes in the direction of rotation by selecting a time with the slider.

Rotation duration:

Value: 10 to 120 s.

2.1.5. Random direction

**Random direction:**

Sets the tunnel orientation randomly for a given time.

Make sure to select at least two directions **by clicking on the corresponding arrows** (to deactivate a direction click on the arrow again) to see the effects of this setting.

This setting can be modified during the session (see Shortcuts).


Keep direction for:

Value: 5 to 60 s.

3. Module

3.1. Session

Once the presets have been selected, click on "**Start the session**" in the bottom right of the screen.



During the session, the user can modify the parameters from the **left side of the screen**. They are not visible to the patient.



In the bottom right of the screen, the user can observe the results of the session live.

3.2. Shortcuts

During the session, the shortcut list is found by clicking on the Xbox controller icon in the upper right corner of the screen.

3.3. Results

Once the session is over, you can view the results.





3.3.1. Summarized results

Summarized results are not available for this module.

3.3.1. Graphs and report

To find the detailed results and to generate a report, click on the histogram icon.



Session details
25/04/2024 11:04



Optical Flow (Advanced)



Parameters

Results

Notes

Name	Initial value	Final value
Environment choice	Advanced	Advanced
Speed	15.00 m/s	15.00 m/s
Automatic inversion	✓	✓
Forward flow duration	10.00 s	10.00 s
Backward flow duration	8.00 s	8.00 s
Random direction	✓	✓
Keep direction for	15.30 s	15.30 s
Opacity filter	0.00	0.00
Countdown	✗	✗
Duration planned	0.00 s	0.00 s
Furrows	✓	✓
Rotation	0.00 °/s	0.00 °/s
Reverse rotation	✓	✓
Rotation duration	90.00 s	90.00 s
Enable fixed point	✗	✗
Spot size	30.00 °	30.00 °

 Start session
with initial values

 Start session
with final values