



Virtualis

Instructions for Use

StaticVR

Table of contents

1.	Introduction.....	4
1.1.	About this IFU.....	4
1.1.1.	Manufacturer	4
1.1.2.	Label	5
1.2.	Product description	7
1.3.	Principles of operation of the device	7
1.4.	Intended use.....	9
1.5.	Clinical benefits	9
1.6.	Indications for use	9
1.7.	Contraindications	9
1.8.	Intended user	10
1.9.	Caution symbols	10
1.10.	Warnings and precautions.....	11
1.11.	Side effects and residual risks	13
2.	General technical specifications.....	14
2.1.	Minimum technical requirements.....	14
2.2.	StaticVR key functional elements.....	15
2.3.	Module compatible accessories	16
3.	Operating instructions.....	19
3.1.	System start / stop	19
3.1.1.	Start	19
3.1.2.	Plugging in the static rehabilitation force plates (old & new generation)	21
3.1.3.	Stop.....	23
3.2.	Before starting a session	24
3.3.	Patient setup	25
3.3.1.	General patient setup.....	25
3.3.2.	Holding the controllers.....	27
3.4.	Starting a module	29
3.5.	Session settings	31
3.6.	Force plate tare and statokinesigram	32
3.7.	During a session.....	34
3.8.	Accessing the results	35
3.9.	Daily maintenance.....	40
4.	Care and maintenance	41

- 4.1. Cleaning procedure 41
 - 4.1.1. Computer cart 42
 - 4.1.2. Computer tower 42
 - 4.1.3. Other components 42
- 4.2. Caring for the headset..... 45
- 4.3. Maintenance requirements..... 45
- 4.4. Decommissioning and disposal 46
- 5. Annexes 47
 - 5.1. Annex 1: Modules included in the StaticVR package 47
 - 5.2. Annex 2: StaticVR modules and their accessories..... 51

1. Introduction

The system is fully installed and checked on site by authorized technicians.

Installation includes:

- System assembly
- Verification of main functionalities

1.1. About this IFU


Non-contractual photos.

Photos shown are for illustrative purposes only. Differences between photos and reality may occur as a result of updates or changes in equipment without impact on safety nor performances.

The abbreviation VR is used for Virtual Reality (VR).

1.1.1. Manufacturer

Virtualis is the legal manufacturer of the medical device software:

<p>Manufactured by:</p> <p>VIRTUALIS</p> <p>Virtualis 78 Allée John Napier Immeuble ATRIUM, 34000 Montpellier France ☎ +33 9 80 80 92 91 ✉: contact@virtualisvr.com</p>	<p>CE marking's date of obtention: 2018</p> <p></p>
---	---

Virtualis is responsible for placing this system, composed of software and various devices, on the European market:


Virtualis
78 Allée John Napier Immeuble ATRIUM,
34000 Montpellier
France
☎ +33 9 80 80 92 91
✉: contact@virtualisvr.com

1.1.2. Label

1.1.2.1. System label


The system is identified with a label affixed to the computer.

Depending on the manufacturing site, the system can be labeled with a label following these templates:



The system can be delivered in several configurations table:

- Desktop
- Laptop
- Wired headset
- Static rehabilitation force plates (old generation, new generation)



1.1.2.2. Software label

The software is labeled as a medical device with the following label displayed on the Patient manager’s home page:



1.1.2.3. Label symbols

Symbol	Description
	Manufacturer of the device
	Manufacturing date
	Part number
	Model number
	Serial number
	Consult electronic instructions for use
	Reference to the warnings and/or precautions to take described in this document
	Unique Device Identifier
	Medical Device
	Humidity limitation
	Atmospheric pressure limitation
	Temperature limit
	The CE-mark indicates that the Virtualis product meets the requirements of the Medical Device Directive 93/42/CE.

1.2. Product description

StaticVR consists of a software package installed on a computer paired with the necessary virtual reality accessories and static rehabilitation force plates with compatible foams.

There are 14 modules included in the StaticVR range. The range includes both assessment and rehabilitation software.

The modules and their availability can be found in Annex 1.

The software is intended to be used for patient assessment and rehabilitation.

1.3. Principles of operation of the device

The solution developed by Virtualis is based on virtual reality technology, the principle of which is to visualize images in a 3D environment using a headset, and to use positioning sensors.

Virtualis develops software to create a virtual environment conducive to patient rehabilitation, thanks to the virtual reality products available on the market.

The software's functionalities (called modules) provide to the practitioner exercises that he/she chooses and configures according to the rehabilitation program he/ she needs for the patient.


The Virtualis solution contains also Patient manager software enables practitioners to save patient data, to use Virtualis modules through parameterized sessions, and to store the results of these sessions. Once saved, these sessions can be consulted or displayed in printable reports.

The practitioner starts a module/exercise via the Patient Manager and set the parameters according to his/her intended rehabilitation program. These inputs parameters are sent to the Module. Use of force plates require calibration before use.

Once the practitioner launches the session/exercise, the virtual reality images are displayed in the headset thanks to the virtual reality software (SteamVR) provided with the virtual reality accessories. Then the patient can see the image displayed in the headset. This image is also available for the practitioner on the screen.

During the session/exercise, the module will get the posturographic data of the patient in real time from force measurement force plates. The virtual reality accessories send data to the virtual reality software (SteamVR) in accordance with their intended use. This data is sent to the module which uses them to provide output parameters and calculated results.

At the end of the session, the output parameters and the results are sent to the Patient Manager.



In addition, remote controls are available by the practitioner to adapt the parameters of the exercise/module while remaining at the patient side during the session.

1.4. Intended use

The device provides quantified assessments of rehabilitation and exercises to enable the healthcare professional to adapt the rehabilitation protocol to the patient's needs.

1.5. Clinical benefits

1. Assess patient disorder by using VR and/or static rehabilitation force plates.
2. Rehabilitate patient function by using VR and/or static rehabilitation force plates.

1.6. Indications for use

Assessment and rehabilitation, including:

- Balance and vestibular,
- Neurology,
- Musculoskeletal.

1.7. Contraindications

Contraindications for the use of virtual reality are not limited to, but may include:

- Patients with photosensitive epilepsy
- Pregnancy
- Per HTC Vive hardware manufacturer: young children are not recommended to use the product, however, there is no age limit recommendation.
- Patients with long-lasting cyber malaise or motion sickness.

The use of virtual reality for therapeutic purposes is not recommended for patients who have not undergone a prior vestibular assessment.

1.8. Intended user

Any Healthcare Professionals involved in rehabilitation/assessments including but not limited to physiotherapists, occupational therapists, speech and language therapists and audiologists, orthotists and prosthetists, clinical psychologists, physical medicine and rehabilitation doctors, and rehabilitation nurses.

The target patient population includes individuals of 15 years and above, excluding epileptic patients and pregnant women.

1.9. Caution symbols



WARNING

This symbol precedes a warning about conditions or practices that may present danger to the patient and/or user.



CAUTION

This symbol precedes a warning about conditions or practices that could result in damage to the equipment.



RECOMMENDATION

This symbol indicates advice on the use of the system, with no impact on the health or safety of users, patients or the environment.

1.10. Warnings and precautions

Immersion in Virtual Reality is a powerful tool, especially for stimuli that can induce sensory conflicts.

WARNING



These stimulations can potentially cause certain disorders: vagal discomfort, epileptic seizures, migraines, vomiting, malaise, dizziness, syncope etc.

This type of re-education must be approached progressively, particularly in Virtual Reality where the stimulation is "powerful".

WARNING



In case of a software malfunction, the screen could flash in the headset and potentially cause epileptic seizures: **do not use VR with epileptic patients.**

WARNING



Monitor the patient to prevent falls.

WARNING



It is recommended that the clinicians who administer test to consider the learning effect.

RECOMMENDATION



As postural reactions can be intense, place the patient in a secure environment and stay close to them throughout the session to anticipate any loss of balance or discomfort caused by the use of virtual reality.

RECOMMENDATION



Increase the duration and intensity of stimulation very gradually after an initial short session to ensure the patient's tolerance to this type of stimulation.

The use of virtual reality must respect the neurosensory capacities of the patient being treated.

**RECOMMENDATION**

Virtual reality for therapeutic purposes can be offered as a complement to conventional vestibular rehabilitation. In current available evidence, virtual reality rehabilitation has not shown to be superior to conventional vestibular rehabilitation.

Motion sickness is treated by "habituation", so you need to recreate the symptoms experienced during daily activities.

**RECOMMENDATION**

Monitor patient for cyber malaise or motion sickness symptoms and their reaction to the virtual reality simulation.

**WARNING**

It is essential to stop the session when the first symptoms appear, generally "sweating".

Define a working safe work area of about 4 m² to allow for risk-free movements.

Take a 10-to-15-minutes break every 30 minutes of use.

It's up to the professional to "dose" immersion so as not to provoke neurovegetative symptoms. This type of symptom can intensify in the hour following the session.

According to HTC VIVE, the VR accessories required to use the software may emit radio waves that can interfere with the operation of nearby electronic devices. If you have a pacemaker or other implanted medical device, do not use the product until you have taken advice from your doctor or the manufacturer of your medical device.



Any serious incident should be notified in writing to <https://virtualis-us-dgs.happyfox.com> for the USA or service@interacoustics.com and the Member State in which the user and/or patient is established.

1.11. Side effects and residual risks

**WARNING**

Transient increase in symptoms including nausea, emesis, dizziness, imbalance, motion sickness, musculoskeletal discomfort, symptoms of hypotension.

**WARNING**

Falls or near falls resulting in injury including fractures.

**WARNING**

Some physicians may want to delay exercises during the early postoperative stage because of risk of bleeding or cerebrospinal fluid leak.

**WARNING**

Dizziness and imbalance as side effects of the exercises could increase psychological distress in some individuals.

**WARNING**

If the software has some performance issues, the patient could experience motion sickness as a result.

2. General technical specifications

2.1. Minimum technical requirements

In order to install and use our virtual reality applications, we recommend a configuration equal to or higher than the system requirements:


Technical Minimum Requirements – Computer	
CPU	i7-9700 (3.5GHz) or greater
GPU	Nvidia GTX 2060 or greater
Operating system	Windows 10 or higher
VRAM	6 Go minimum
RAM	16 Go minimum
Storage	512 Go minimum

Technical Minimum Requirements – Static rehabilitation force plates	
Connection	2x USB port (1 per plate)
Weight capacity	200kg maximum
Dimensions (each plate)	34,5 x 45 x 2,2 cm (LxWxH)
Weight	40000g
Sensitivity	20gr
Accuracy (CoP)	+/- 0,55 mm
Frequency	75Hz
Bilateral CoP	Each foot independently

Technical Minimum Requirements – Other	
Screen	Resolution: 1920x1080 or greater
VR headset	HTC Vive series
Workstation	64 x 60 x 170 cm (W x L x H)

2.2. StaticVR key functional elements

The medical device software consists of an application (called “Patient Manager”) which enables the use of software modules that the therapist selects according to his/her assessment or rehabilitation needs. Software modules are marketed in packages as illustrated in the following figure:



The key functional elements of the device are:

- **Patient manager that includes:**
 - Create patient file
 - Set the modules parameters
 - Start a module
 - Provide results
 - Generate reports
 - Create chaining
 - Administrative management

- **Modules:**
 - Display VR images
 - Edit parameters
 - Start the simulation
 - End and exit the exercise
 - Calculate and display results
 - Access patient’s center of pressure data

2.3. Module compatible accessories

Compatible accessories vary depending on the selected module. They are the following:

- VR headset
- Handheld controllers
- Xbox controller
- Secondary screen (optional)
- Static rehabilitation force plates
- Foam frame
- Airex Balance pad Elite foam pad or similar models

**CAUTION**

Do not use any other devices than the ones provided by Virtualis.

The accessories and equipment may vary depending on the chosen package or the options you may have purchased. See Annex 2 for more information on the equipment available and which modules they are compatible with.

The software runs on a computer and screen with adequate performances to support virtual reality.

The interactive virtual reality software modules use head mounted display virtual reality hardware with a virtual reality- ready PC, to deliver the Virtualis VR Software to the patient.

The computer can be installed on a cart to provide the user with a control station, composed of:

- Computer Monitor
- Storage Basket
- Virtual Reality-Ready Computer Tower
- 4-Port USB Hub
- Link Box
- A second screen (optional)
- Static rehabilitation force plates






FIGURE 1 - COMPUTER CART EXAMPLE

CAUTION



Prevent the accessories from falling by setting them in the storage basket, or on a sufficiently large table. For setups without a cart, reserve a space of 90 x 60 cm minimum to store the accessories.

CAUTION



The computer provided is only intended for professional use of the Virtualis software. Do not install unauthorized software or modify the computer's initial configuration.

Remote controllers are supplied to be used by the practitioner as a remote control to navigate software menus from a distance while at the patient's side:

- A wireless Xbox controller:



FIGURE 2 - XBOX CONTROLLER

The software package is intended to operate with compatible virtual reality accessories verified and validated by Virtualis. These virtual reality accessories are intended to be held by the patient while they perform the exercises displayed by the software. The choice of accessories depends on the software module and its intended purpose.

Virtual Reality Headset:

- Headset Strap (adjustable)
- Headset Adjustment Dial
- VR Lenses
- VR Headset and Face Cushions, that are in direct contact with the patient
- Earphones



FIGURE 3 - VIRTUAL REALITY HEADSET EXAMPLES

2 controllers are provided with their charging cables.

Some software modules can be used by the patient without the VR headset, viewing the environment via a 2nd optional screen.

Some software modules are intended to be used with **memory foam blocks** (see Annex 2 for more information). We recommend the use of **Airex Balance Pad elite** foam blocks, or products with similar characteristics.

Some of the software modules are compatible with the static force plates and with another Virtualis device called Dynamic balance board for MotionVR package(see **Annex 2**).

3. Operating instructions

3.1. System start / stop

3.1.1. Start

Follow the steps below:

- 1) Connect the main plug to the main socket.
- 2) Turn ON the computer
- 3) Start the link box in order to connect the headset to the computer or turn ON the router and turn ON the headset (depending on your configuration).



FIGURE 4: EXAMPLE: HTC VIVE PRO AND FOCUS VISION HEADSET'S "POWER" BUTTON


- 4) The system starts up and displays the SteamVR software or Vive Hub on the screen:

Virtual reality accessory state

List of virtual reality accessories



Virtual reality accessory state



CAUTION



If the system is connected to the Internet, SteamVR software can be updated automatically.


If the update affects the Virtualis software operation, please contact the customer service.

USA: <https://virtualis-us-dgs.happyfox.com>


Other countries: service@interacoustics.com.

3.1.2. Plugging in the static rehabilitation force plates (old & new generation)

- 1) Place the left and right plates into the foam frame



- 2) Plug both USB cables into the USB-C Hub.



- 3) Plug the USB-C Hub into the computer's USB-C port.



FIGURE 8: EXAMPLE: INSTALLATION COMPLETE

3.1.3. Stop



CAUTION

To avoid any risk of data loss, always switch off the system and wait for it to come to a complete standstill before disconnecting the power cable.



RECOMMENDATION


Make sure the computer is properly turned off at the end of each day.

Follow the steps below:


- 1) Close/ Exit Virtualis Patient Manager software.



- 2) Close/ Exit SteamVR.



- 3) Turn off the system from Windows menu.



- 4) Turn off the screen.
- 5) Turn off the link box (depending on your configuration).
- 6) Disconnect the main plug to the main socket.

3.2. Before starting a session





RECOMMENDATION

Familiarize yourself with the modules before using them with a patient.




RECOMMENDATION

Check the battery level of the VR accessories and charge them if need be.

- 1) Launch the Virtualis software. 
- 2) Launch the SteamVR application. 
- 3) Check that the headset and the base stations are detected by the SteamVR or Vive Hub application.

The accessories appear in blue if they are properly detected.



EXAMPLE VR ACCESSORIES PROPERLY DETECTED BY STEAMVR

- 4) Select a patient in the Virtualis Patient Manager.
- 5) Check that the patient information displayed matches your patient's. Edit it if necessary.



3.3. Patient setup

3.3.1. General patient setup

**RECOMMENDATION**

The patient setup differs depending on the module you have chosen. Check the module's tutorial to find it.

- 1) Ensure the patient is properly situated and has enough room to move freely without risk (a minimum of 4 m² is required).
- 2) The force plates' capacity is 200kg maximum.
- 3) Always stay nearby to prevent them from losing their balance or falling.
- 4) Help the patient put on the VR headset. It should be tightly but comfortably fixed on the patient's head. It can be worn with glasses.
- 5) Help the patient equip any required accessory depending on the module you have chosen.
- 6) Show the patient how to hold the controller(s) to ensure they will not fall.
- 7) Give them all the necessary indications prior to starting the session:
 - The module you have chosen,
 - Its objective(s),
 - The movements they will have to perform.
- 8) Ensure the patient keeps their hands off the headset's sensors.
- 9) Put the Static force plates into the foam frame to keep them into place and keep them from slipping.
- 10) Help the patient onto the force plate.

**RECOMMENDATION**

The medial malleolus of each of the patient's feet must be centered directly on the horizontal line of the force plate.

- 11) Use an Airex Balance pad Elite or similar models cushion if the module requires this option: place the cushion on the force plates, aligning its markings with those of the foam frame:





FIGURE 9: EXAMPLE STATIC FORCE PLATES & FOAM FRAME & BALANCE PAD

Positioning the patient

The device must be correctly turned on and initiated and computer software operational prior to patient positioning on the device.

When standard foot placement is required, position the patient either directly on the force plates or the Airex cushion as follows:

1. Help the patient step onto the force place or Airex cushion, facing forward
2. Center the patient's feet on the force plate or Airex cushion. The medial malleolus of each foot should be centered directly over the center horizontal (wide) line on the force plates or aligned with the horizontal lines on the Airex cushion.



Cushion and foam



Feet placement foam


3.3.2. Holding the controllers

Show the patient how to hold the controller(s) to ensure they will not fall.

- Vive Pro controllers:



System button

Make sure the patient keeps their fingers off the system button , which can quit a module.



- Vive Focus series controllers



Menu button




VIVE button

Make sure the patient keeps their fingers off the Menu button  and the VIVE button , which can quit a module.

3.4. Starting a module

Refer to the Patient Manager’s User manual for more information.




RECOMMENDATION

If you are unable to log in, contact the administrator to unlock the account and/or change the password.


Once the patient is selected in the Patient manager, properly set up in the room and the VR accessories are properly detected:

- 1) Click on **“New session”** or **“Start a session”** in the Virtualis software.




- 2) Choose a module in the list.


- Use the **“Search”** bar to find a specific module.




- Click on “**Module choice**” to filter assessment or rehabilitation modules.



- You can also filter the modules by pathology.




- 3) Read the module’s description by clicking on the (i) icon and check out the tutorials (video and document). The written documentation can be found in a dedicated folder on the computer.



THE MODULE VISUAL MOTION SENSITIVITY IS USED AS AN EXAMPLE

- 4) Quit the module description

5) Click on "Launch".



3.5. Session settings

Once a module has been selected:

- You can adapt the session to your patient's needs and to your objectives by adjusting the settings.
- The number of parameters available varies depending on the chosen module.






FIGURE 5 : EXAMPLE OF THE SETTINGS AVAILABLE FOR THE MODULE "WEIGHT BEARING"

3.6. Force plate tare and statokinesigram

To use the softwares with the static force plates, follow these steps:

Note : The display might vary depending on the module, but the steps remain the same regardless.

1. Click on the "Tare" button before letting the patient step up onto the force plate.




2. Help the patient onto the force plate.
3. If necessary click on the arrow in the upper right corner of the screen to display the force plates' settings.



4. If necessary modify the force plates' settings at any point during the session
5. Modify the force plates' settings at any point during the session

3.7. During a session



A list of **shortcuts** can be found in most of the modules by clicking on the keyboard or Xbox controller



icon **in the top right corner of the screen.** You can scroll through different shortcuts

options by clicking on the arrows:

In most modules, VR environment settings are available and can be adjusted from the **left side of the screen.**

Static force plates settings can be adjusted from the **right side of the screen.**

You can have access to remote controls using the Xbox controller, the footswitch, the controllers, etc.



CAUTION

The patient could leave the session by themselves by pressing the wrong button on the controllers. **Make sure the patient is holding the accessories properly.**

The module pauses and a window offers to confirm or cancel the requested action. The patient must cancel the action by clicking on “cancel” to return to the virtual reality environment.

3.8. Accessing the results

At the end of a session:




RECOMMENDATION

Not all modules have results to access (summarized or detailed ones).

- **Accessing the summarized results:**

1) Access the session's summarized results by clicking on "**History**" in the Patient Manager.



2) Select the session whose results you wish to see.

3) You access a page showing the initial and final parameters.

Session details
27/05/2025 15:27
Weight Bearing

Parameters
Results
Notes

Name	Initial value	Final value
Starting foot	Left	Left
Max left foot pressure	70 %	70 %
Max right foot pressure	70 %	70 %
Tolerance level	20 %	13 %
Hold time	3.00 s	3.00 s
Session time limited	×	×
Reward step	10	10

Start session
with initial values

Start session
with final values


4) Click on “Results”. (Some modules do not have summarized results).

Session details
27/05/2025 15:27
Weight Bearing


Parameters
Results
Notes

Name	Result
Hold time succeeded	12
Hold time attempt	14
Hold time score	86.00 %
Elapsed time	61.00 s

5) Click on “Notes” if you want to keep Notes




- **Accessing the detailed results and the report:**

1) Click on the graph icon .



2) Select a type of graph.



3) Click on “**Chart options**”  to choose which information to display on the graph.

More options

Include properties

Unselect all

Results

Hold time succeeded	+
Hold time attempt	+
Hold time score	% +
Elapsed time	s +

Parameters

Starting foot	-
Max left foot pressure	% -
Max right foot pressure	% -
Tolerance level	% -
Hold time	s -
Session duration	s -
Reward step	-

4) Click on “Add to report” to add the graph to a report.

Report
Edit and visualize reports for Demo Physio Event

Edition PDF

Chart options

Weight Bearing (27/05/2025)

Legend:
Hold time succeeded (purple)
Hold time attempt (blue)
Hold time score (%) (green)
Elapsed time (s) (pink)

Metric	Value
Hold time succeeded	12
Hold time attempt	14
Hold time score (%)	86
Elapsed time (s)	61

+ Add to report

5) Click on “PDF” to see the generated document.

3.9. Daily maintenance

- Clean the headset and accessories daily, following the instructions in the product manufacturer guide.
- After each session with a patient, it should be ensured that there is no contamination on the equipment and accessories. General precautions must be observed to avoid the transmission of infections and diseases between patients. **Sanitize the accessories between each patient, including the headset’s lenses.** (See 4.1 Cleaning procedure for more information).
- Perform a check each morning and evening to make sure the equipment is clean, working properly and fully charged.


Icon	Status	Description
	Low battery	The accessory’s battery is low.

TABLEAU 1: LOW BATTERY ICON IN STEAMVR

- Charge all the accessories at the end of each day.

4. Care and maintenance



WARNING

To prevent data loss or theft, choose a secure Patient Manager password and make sure to keep it private.

See the Patient Manager's manual for more information.



RECOMMENDATION

Perform software updates as soon as they are available.

4.1. Cleaning procedure



CAUTION

The following points are to be used as general guidelines; **always refer to product manufacturer guide for cleaning instructions for specific products or components.**



WARNING

Always use protective gloves when cleaning and disinfecting equipment.



RECOMMENDATION

Dust the station regularly, using a soft, damp antistatic cloth soaked in cleaning solution.

4.1.1. Computer cart

**WARNING**

Switch off and disconnect the system from the main power supply before cleaning it.

4.1.2. Computer tower

**WARNING**

Do not spray or pour liquid cleaners on the central unit. Liquid may penetrate the ventilation holes. Use only a soft cloth dampened with cleaning solution.

The central unit can be cleaned with a soft cloth soaked in cleaning and disinfectant products commonly used in hospitals, including isopropyl alcohol and neutral cleaning solutions.

Substances not to be used:

1. thinner, benzene, alcohol
2. insecticides or other volatile elements
3. organic solvents such as acetone and toluene

4.1.3. Other components

**WARNING**

Avoid splashing liquids on electrical components.

Do not use abrasive products.

Use isopropyl alcohol or neutral cleaning solutions.

4.1.3.1. Screen

**WARNING**

Do not spray or pour liquid cleaners onto the screen vents.

**WARNING**

When cleaning the screen, take care not to spill any liquid on the central unit.

Lightly wipe dirt from the screen surface with a soft cloth soaked in a neutral cleaning solution.

Substances not to be used:

- Thinner, benzine, alcohol
- Insecticides or other volatile elements
- Organic solvents such as acetone and toluene
- Avoid long-term contact with rubber or vinyl products.

Do not rub the surface of the LCD panel with a rough object, and do not press hard on the surface of the LCD panel.

4.1.3.2. Virtual reality accessories

**RECOMMENDATION**

Please refer to the instructions provided with the accessories.

Headset cushion:

Clean with a nonabrasive antibacterial/disinfecting wipe between uses *or* per your facility and institutional policy for cleaning non-porous foam. It is also possible to use an additional protective barrier between the patient and VR headset components.

4.1.3.2.1. Headset**RECOMMENDATION**

Clean the headset between each patient.

a) Cleaning the face cushion (VR cover)

See: https://www.vive.com/us/support/vive-pro2/category_howto/cleaning-the-face-cushion.html and https://www.vive.com/us/support/vive-pro2/category_howto/replacing-the-face-cushion.html#replacing-the-face-cushion for more information.

- 1) To remove the VIVE Pro II's face cushion, peel off starting from both ends of the face cushion until it detaches from the headset.
- 2) When cleaning the face cushion, remember to:
 - Dampen a clean smooth cloth with cold water and gently wipe the area you want to clean.
 - Air dry it at room temperature. Don't tumble dry, iron, or expose it to direct sunlight.
 - Don't scrub, wring, or bleach it.
 - Don't dip or soak it in water.
- 3) To replace the face cushion, insert the tab on the face cushion into the slot on the headset.
- 4) Align the hook and loop strips.

b) Cleaning the exterior of the headset and the controllers

See: https://www.vive.com/us/support/vive-pro2/category_howto/caring-for-your-vive.html for more information.

Wipe them down using a clean microfiber cloth.

c) Cleaning the headset lenses

See: https://www.vive.com/us/support/vive-pro2/category_howto/cleaning-the-headset-lens.html for more information.

- 1) Make sure the cleaning cloth provided in the box is dry.
- 2) Wipe the lenses in a circular motion from the center to the outer edges of the lenses.

**CAUTION**

Do not scratch the lenses or disassemble other parts of the headset.

**RECOMMENDATION**

You can use an air pump to blow away dust and particles from the headset lenses.

4.1.3.2.2. Other accessories

Refer to the manufacturer's care instructions for more information.

4.2. Caring for the headset

Refer to the VIVE care instructions:

https://www.vive.com/us/support/vive-pro2/category_howto/caring-for-your-vive.html

4.3. Maintenance requirements

For any maintenance inquiry, such as log files to be checked or possibly cleared, contact authorized personnel.

4.4. Decommissioning and disposal

If decommissioning or disposal of the software is necessary, the medical data should be exported as PDF files from the Patient Manager.

Then, the software should be uninstalled, and the local data folder should be removed from the computer.

To ensure that the data is completely deleted, it is recommended to format the hard drive.

5. Annexes

5.1. Annex 1: Modules included in the StaticVR package

Table legend:
I: Included
X: Excluded

Modules	Assessment / Rehabilitation	Clinical application
Dynamic Analysis	Assessment	This module quantifies weight distribution between the right and left sides.
BirdVR	Rehabilitation	Pathologies combining the need to modulate support and balance. Dissociation of girdles . Balance disorders in general. Double-task executive disorders.

Modules	Assessment / Rehabilitation	Clinical application
CTSIB VR	Assessment	Evaluation of sensory interaction in balancing reactions.
Supermarket Scrolling	Rehabilitation	Balance disorders. Vestibular pathologies.
Unilateral Stance	Assessment	Balance problems. Risk of falling.
MazeVR	Rehabilitation	Impaired postural adaptation and loss of stability. Cognitive, attentional, concentration and problem-solving disorders.
LOS	Assessment	Measure of a person’s limits of stability by giving an angle measurement of the tilt cone.

Modules	Assessment / Rehabilitation	Clinical application
LOS Rehab	Rehabilitation	Rehabilitation of the patient's balancing strategy according to the limits of the angle of his oscillation cone.
MatchingVR	Rehabilitation	Any musculoskeletal or neurological pathology of the upper limbs and spine. Balance and posture disorders. Memorization.
Pursuit CoP	Rehabilitation	Functional rehabilitation of body weight transfer and pressure in orthopedic, neurological, rheumatological or traumatic contexts of the lower limb.
Weight Distribution	Assessment	Posturographic assessment of plantar pressure distribution.

Modules	Assessment / Rehabilitation	Clinical application
SkiVR	Rehabilitation	Orthopedic, neurological, rheumatic functional re-education of the lower limb. Center of gravity displacement work.
Weight Bearing	Rehabilitation	Rehabilitation of a patient's pressure distribution and balancing strategy according to body weight distribution.
Visual Motion Sensitivity	Assessment	Vestibular pathologies, balance disorders, vertigo.

5.2. Annex 2: StaticVR modules and their accessories

Table legend:
I: Included
X: Excluded

Modules	Patient required VR Accessories									Partient optional VR Accessories	Other accessories
	VR Headset	VR Headset and Controller(s)	VR Headset and Tracker(s)	VR Headset and Controller(s) and Tracker(s)	VR Headset and Xbox controller	VR Headset and keyboard	VR Headset or 2nd screen and Force plate	StaticVR Force plate	MotionVR dynamic platform		
Dynamic Analysis	X	X	X	X	X	X	X	I	A	X	X
BirdVR	X	X	X	X	X	X	I	I	A	Controller	Second screen
CTSIB VR	X	X	X	X	X	X	I	I	A	X	Airex Balance Pad elite memory foam
Supermarket Scrolling	X	X	X	X	X	X	I	I	A	X	Second screen

Modules	Patient required VR Accessories									Partient optional VR Accessories	Other accessories
	VR Headset	VR Headset and Controller(s)	VR Headset and Tracker(s)	VR Headset and Controller(s) and Tracker(s)	VR Headset and Xbox controller	VR Headset and keyboard	VR Headset or 2nd screen and Force plate	StaticVR Force plate	MotionVR dynamic platform		
Unilateral Stance	X	X	X	X	X	X	I	I	A	X	X
MazeVR	X	X	X	X	X	X	I	I	A	X	Second screen
LOS	X	X	X	X	X	X	I	I	A	X	X
LOS Rehab	X	X	X	X	X	X	I	I	A	X	Second screen

Modules	Patient required VR Accessories									Partient optional VR Accessories	Other accessories
	VR Headset	VR Headset and Controller(s)	VR Headset and Tracker(s)	VR Headset and Controller(s) and Tracker(s)	VR Headset and Xbox controller	VR Headset and keyboard	VR Headset or 2nd screen and Force plate	StaticVR Force plate	MotionVR dynamic platform		
MatchingVR	X	A	A	X	X	X	A	I	A	X	Second screen
Pursuit CoP	X	X	X	X	X	X	I	I	A	X	X
Weight Distribution	X	X	X	X	X	X	X	I	A	X	X
SkiVR	X	X	X	X	X	X	I	I	A	X	Second screen

Modules	Patient required VR Accessories									Partient optional VR Accessories	Other accessories
	VR Headset	VR Headset and Controller(s)	VR Headset and Tracker(s)	VR Headset and Controller(s) and Tracker(s)	VR Headset and Xbox controller	VR Headset and keyboard	VR Headset or 2nd screen and Force plate	StaticVR Force plate	MotionVR dynamic platform		
Weight Bearing	X	X	X	X	X	X	I	I	A	X	Second screen
Visual Motion Sensitivity	X	X	X	X	X	X	I	I	A	X	X