

# PA5

## Feature Overview



**Interacoustics<sup>®</sup>**

leading diagnostic solutions

### **Conditioned Orientation Response**

Two PA5s can be used for the COR test by using one in each hand. The tester can use either unit to evoke a head turn with the stimulus and then provide visual reinforcement. The single-hand controls on the PA5 allow easy control over frequency and intensity selections for faster testing.

### **Extended Battery Life**

The PA5 is powered by three standard AA batteries. It will turn off automatically after a period of inactivity and reactivate when a stimulation switch is touched.

### **Silent Presentation Buttons**

The two stimuli, light and sound, are controlled silent buttons. Both stimuli can be presented simultaneously by letting one finger touch both buttons. This gives full freedom operating COR, VRA and APR tests.

### **Improved Sensitivity**

The PA5 includes white noise as well as warble and narrow band stimuli. White Noise has been shown to be especially useful for evoking responses in the early months.

### **Audiometric Headphone**

An optional TDH39 headphone may be connected to the PA5 to test more cooperative children who are capable of providing responses. This allows for more accurate hearing threshold assessments. Once connected, the PA5 recognizes the TDH 39 phones and adjusts its calibration accordingly.