

# RhinoStream Tests



**Interacoustics<sup>®</sup>**

leading diagnostic solutions

### **Start Measuring**

Once the RhinoStream hardware and software have been started, and the probe initialized you are ready to start measuring a patient. You have a choice between either using the standard clinical mask including the manometer probe, or screening patients by attaching the probe to the nose adapters. The patient should be instructed in how to breathe correctly when consecutive measurements of both the left and right nose are made. It only takes a short time to get data for your diagnosis if you use the screening procedure.

### **Display of the results according to**

- The International Standardization Committee's Recommendation
- The Japanese Standard