

Science **made** smarter

AD528

Diagnostics
made easy

Easy-to-use
air, bone and speech
audiometer with
comfortable
ergonomics



Interacoustics

Audiometry
~~~~~

Tympanometry

ABR

OAE

Hearing Aid Fitting

Balance



# Diagnostics made easy

## Compact and easy to use

AD528 is the perfect audiometer for you if you are looking for diagnostic power that can be scaled according to your specific test needs. We proudly call it 'Diagnostics made easy':

- Starts in a few seconds. Makes it ideal for the busy clinic or traveling clinician with many power ups during a day
- Easy to use. The simple button layout ensures a short learning curve
- Comfortable ergonomics. With an angled keyboard and convenient keyboard layout, you can comfortably test your patients all day long
- Fast masking. You can quickly apply or cease masking with the keyboard's right wheel

## Flexible test battery

AD528 offers a flexible test battery that includes the Stenger and Modified Hughson-Westlake tests. You can tailor AD528 to your needs by choosing between the additional tests.

## Comfortable ergonomics

Ergonomics is a key aspect of AD528. The angle of the keyboard together with the conveniently positioned tone switch and wheels secure a relaxed working position. This will save your wrists from strain and make it more comfortable to perform measurements all day.


## Have it in your language

You can change the display language to Chinese, Czech, English, Finnish, French, German, Greek, Italian, Japanese, Korean, Norwegian, Polish, Portuguese, Russian, Spanish or Turkish.

You can also choose the keyboard text to be in one of these languages: Brazilian Portuguese, Chinese, English, French, German, Italian, Japanese, Korean, Polish, Russian, Spanish or Turkish.

## Design your own reports and printouts

The optional software suite includes a print wizard for customizable reports and printouts. You can also print directly to a small 3" printer.



Merge and configure input from different sources in the print wizard



**Ideal for traveling**

The small footprint and low weight of 1.5 kg combined with the optional carrying bag and fast start up time make AD528 the natural companion for the clinician on the road.

**Integrate with numerous systems**

You can purchase a sync license for integration with Noah, local EHR systems or OtoAccess® Database.



**50 years of evolution**

Interacoustics has been developing audiometers for 50 years. This is your guarantee for a quality solution where all the benefits of previous audiometers have been compiled to match the needs of a busy clinic that requires reliability every single day.

| Test battery                                                                     | AD528    |
|----------------------------------------------------------------------------------|----------|
| Air, bone and free-field                                                         | x        |
| Speech (CD, MP3, Microphone or Wave)                                             | x        |
| Stenger                                                                          | x        |
| Modified Hughson-Westlake                                                        | x        |
| ABLB                                                                             | Optional |
| SISI                                                                             | Optional |
| Langenbeck                                                                       | Optional |
| Békésy                                                                           | Optional |
| Weber                                                                            | Optional |
| Speech in noise and Binaural speech                                              | Optional |
| QuickSIN                                                                         | Optional |
| Sync license for integration with Noah, local EHR systems or OtoAccess® Database | Optional |

**AD528 is based on 50 years of experience in developing audiometers**



# Science made smarter

## Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

[Interacoustics.com](http://Interacoustics.com)

Interacoustics A/S


Audiometer Allé 1  
5500 Middelfart  
Denmark

+45 6371 3555  
[info@interacoustics.com](mailto:info@interacoustics.com)

[interacoustics.com](http://interacoustics.com)

Go online to  
explore our  
full product  
range


## Related products



**AT235**  
Middle ear analyzer



**AD629**  
Diagnostic audiometer



**Viot™**  
Video otoscope

### Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.



**Interacoustics**