




# Interacoustics Software Suites



---

Increase patient  
satisfaction and clinic  
efficiency



# The Interacoustics suite concept

Interacoustics provides one of the most comprehensive software based solutions for audiologic diagnostics and hearing aid fitting. The Interacoustics Suites integrate various measurement modules in a common yet intuitive and user friendly way to enable viewing, storing and printing of shared data.

## Diagnostic Suite

Supports all standalone, non-PC based audiometers and tympanometers. OtoAccess™ supports all Interacoustics products. This means that all Interacoustics products are compatible with the HL7.

## Titan Suite

Supports the Titan Middle Ear Analyzer

## Affinity<sup>2.0</sup>/Equinox<sup>2.0</sup> Suite:

Supports PC based Audiometry, Hearing Instrument Test and Real Ear with Visible Speech Mapping.

## Additional modules

Simulators, speech from harddrive, TEN(HL) test, high frequency, etc.

EMR  
(Electronic  
medical records)

Combined  
Reporting  
Print Wizard



Titan Suite  
Tympanometry, OAE & Automated ABR



VIOT™ Video Otoscope Suite



Diagnostic Suite  
Supports all Interacoustics stand alone audiometers and tympanometers




Affinity<sup>2.0</sup> / Equinox<sup>2.0</sup> / Callisto™ Suite  
Audiometry, REM, HIT, and Visible Speech Mapping






- Seamless integration
- Easy workflow
- Intuitive interfaces
- Future safe



Merge and configure input from different sources in the Print Wizard



Printout via Diagnostic Suite software and the Print Wizard.

## Combined printing and reporting

- Retrieve patient demographics from OtoAccess™ and NOAH
- Merge audiometric and impedance test data on a single page
- Text field for session observations and conclusions
- Insert bitmaps such as clinic logo
- Individual session dates
- Highly configurable symbols for quick reference

## Unique Print Wizard for combined reports

Interacoustics Software Suites include a common embedded print wizard. This is a flexible tool that supports individual and clinic reporting requirements to streamline clinic workflow. The Print Wizard combines audiometric, impedance and fitting data for full patient reporting. Unlimited customized templates permit useful alternative templates for different test applications. Printing the reports to PDF allows useful integration into most EMR programs.

## EMR and 3rd party integration

Interacoustics Suites are prepared for the demands of modern Electronic Medical Records and 3rd party systems by using well-defined and documented interfaces with standard protocols like XML and GDT formats. New methods will continually be implemented as new trends arise.

## Data storage

The Interacoustics Suites are designed to integrate data through NOAH or our own proprietary OtoAccess™ database. In addition, all instruments may be run in 'stand-alone' mode and then integrating data via PDF to EMR.

## Future safe platform

Interacoustics Suites are future proof platforms serving as the core for all future software based products. Users will benefit from this concept for years to come.

## Supported instruments and modules

Suite	Modules	Instruments		
Diagnostic	Audiometry Tympometry	AS608/ AS608e AD226 AD229b/e AD629 AC33 AC40 MT10 AT235/AT235h AA222/AA222h		
		Titan	Tympometry	Titan
		Affinity <sup>2.0</sup>	Audiometry REM HIT	Affinity <sup>2.0</sup>
		Equinox <sup>2.0</sup>	Audiometry	Equinox <sup>2.0</sup>

Additional optional modules and simulators for the Affinity<sup>2.0</sup> and Equinox<sup>2.0</sup> Suites: SISI, Speech from hard drive, High-frequency audiometry (up to 20kHz), Master hearing aid, Hearing loss simulation built, Multi-frequency and TEN(HL) test. Diagnostic Suite is Windows® 7 and XP compatible.



---

## Interacoustics

Interacoustics is a world leading provider of diagnostic solutions in the field of hearing and balance assessment. Since 1967 we have designed and manufactured our innovative diagnostic solutions for the world of audiology with a constant focus on providing our customers with quality, dependable products.

---

## With you at all times

We operate in over 100 countries worldwide through a carefully selected network of distributors and service centres. Purchasing an Interacoustics product guarantees not only a quality product, but also direct access to our highly professional training and support service.

Read more at  
[interacoustics.com](http://interacoustics.com)

Go online to  
explore our  
full product  
range



The OtoAccess™ database

---

## Interacoustics A/S

Audiometer Allé 1  
5500 Middelfart  
Denmark

T +45 6371 3555

F +45 6371 3522

[info@interacoustics.com](mailto:info@interacoustics.com)

[interacoustics.com](http://interacoustics.com)

