

Interacoustics - Windows Compatibility Overview

Interacoustics Software	Version	Windows XP	Windows 7 (32 bit)	Windows 7 (64 bit)
Affinity Suite	2.0.4 SP2	YES	YES	YES
Titan Suite	1.0.3	YES	YES	YES
Diagnostic Suite	1.0.1	YES	YES	YES
VNG	7.04	YES	YES	YES
Video Frenzel	1.0	YES	YES	YES
TEOAE25	3.04.2	YES	YES	YES
ABRIS	1.06.1	YES	YES	YES
DPOAE 20	1.03.1	YES	YES	YES
ASSR	1.02.1	YES	YES	YES
EP15/EP25	4.01	YES	YES	YES
OtoAccess	1.2.1	YES	YES	YES

Last updated 20-12-2010 - Copyright © Interacoustics A/S

