

# TBS10 Test Box

Validation made portable

The TBS10 test box offers a complete hearing instrument testing solution when used with our HIT software module and the Callisto™, Equinox<sup>2.0</sup> or Affinity<sup>2.0</sup>.

- Small footprint
- Robust performance
- Coupler based hearing aid fitting (RECD inside the TBS10 test box)
- Frequency range 100-10,000Hz
- Built-in test sequences based on IEC and ANSI standards
- Test sessions can be saved and recalled for later direct comparison

Use TBS10 with the Affinity<sup>2.0</sup> or Equinox<sup>2.0</sup> and get a compact & flexible test box solution offering:

- Supported by a battery drain meter
- Telecoil measurements
- RECD measurements with our REM module
- Use of SKS10 Skull Simulator with TBS10

Hearing aid diagnostics and RECD (real-ear-to-coupler difference) measurements in a compact modular solution.



The Callisto™

The TBS10 Test Box



The Equinox<sup>2.0</sup>



The Affinity<sup>2.0</sup>



# Science made smarter

## **Interacoustics is more than state-of-the-art solutions**

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

[Interacoustics.com](http://Interacoustics.com)

**Interacoustics A/S**

Audiometer Allé 1  
5500 Middelfart  
Denmark

+45 6371 3555  
[info@interacoustics.com](mailto:info@interacoustics.com)

[interacoustics.com](http://interacoustics.com)

Go online to  
explore our  
full product  
range



**Interacoustics**